



## LMA<sup>®</sup> Gastro Airway with Cuff Pilot Technology

### Proactive airway management for endoscopy procedures

#### Separate gastric and airway access


Hypoxemia due to respiratory depression or airway obstruction is a known risk associated with endoscopic procedures, with studies showing that hypoxemia can occur in 11–50% of cases.<sup>1-3</sup> The LMA<sup>®</sup> Gastro Airway with Cuff Pilot Technology from Teleflex is the first laryngeal mask specifically designed to enable clinicians to proactively manage their patients' airways while facilitating direct endoscopic access via the integrated endoscope channel. With the airway in place, clinicians can monitor end-tidal CO<sub>2</sub> for patient safety.

#### Designed to support patient safety

The single-use LMA<sup>®</sup> Gastro Airway is designed for patient comfort with a silicone cuff that is soft and flexible, and conforms to the patient anatomy to create an effective oropharyngeal seal.<sup>4,5</sup> Silicone cuffs have been shown to reduce risk of sore throat<sup>6</sup> and may achieve higher seal pressures<sup>4</sup> compared with PVC cuffs.

#### Integrated Cuff Pilot Technology

Cuff Pilot Technology is an integrated, cuff pressure indicator that constantly monitors cuff pressure, detecting



changes at a glance resulting from fluctuations in temperature, nitrous oxide levels, and

movements within the airway. Cuff Pilot Technology was developed to support clinicians in avoiding the known risks of cuff hyperinflation, which include sore throat, dysphagia, an increased risk of aspiration due to leakage around the cuff, and hypoglossal, lingual or recurrent laryngeal nerve palsies.<sup>7-10</sup>

#### Benefits



##### Clinician

Inspires confidence by supporting end-tidal CO<sub>2</sub> monitoring and direct endoscopic access



##### Institution

Designed to help reduce airway-related complications during endoscopy procedures



##### Patient

Silicone cuff and integrated cuff pressure monitoring for patient comfort

# LMA® Gastro Airway with Cuff Pilot Technology



### Endoscope channel


Enables an endoscope to be passed through the device under vision

### Adjustable holder and strap

Maintains the device in a neutral position during endoscope manipulation

### Integrated bite block

Reduces the potential for damage to, or obstruction of, the airway tube or endoscope due to biting



### Silicone airway tube and cuff

Designed for smooth insertion and patient comfort

### Cuff Pilot Technology

An integrated, single-use cuff pressure indicator that constantly monitors cuff pressure

### Yellow zone

0–40 cm H<sub>2</sub>O

### Green zone

40–60 cm H<sub>2</sub>O

### Clear zone

60–70 cm H<sub>2</sub>O

### Red zone

...70+ cm H<sub>2</sub>O

## LMA® Gastro Airway with Cuff Pilot Technology fast facts

Designed to provide control of a patient's airway while enabling direct access to the esophagus and upper gastrointestinal tract in patients undergoing endoscopic procedures.

In an initial study of 292 patients, the LMA® Gastro Airway was placed and an endoscope inserted successfully in 99% of patients.<sup>11</sup>

Once placed, the LMA® Gastro Airway facilitates end-tidal CO<sub>2</sub> monitoring throughout the procedure to support patient safety.

## LMA® Gastro Airway with Cuff Pilot Technology

ITEM NUMBER	MASK SIZE	PATIENT WEIGHT	INTRACUFF PRESSURE*	MAXIMUM ENDOSCOPE SIZE (OD)
1E5030	3	30–50 kg	60 cm H <sub>2</sub> O	14 mm
1E5040	4	50–70 kg	60 cm H <sub>2</sub> O	14 mm
1E5050	5	70–100 kg	60 cm H <sub>2</sub> O	14 mm

\* Cuff Pilot Technology recommendation for cuff pressure: green zone pressure range = 40 – 60 cm H<sub>2</sub>O for all sizes. Note: The LMA® Gastro Airway does not facilitate intubation.



[teleflex.com/lmagastro](https://teleflex.com/lmagastro)

### References.

1. Qadeer MA, Rocio Lopez A, Dumot JA, Vargo JJ. *Dig Dis Sci.* 2009;54(5):1035-1040.
2. Cote GA, Hovis RM, Ansstas MA, et al. *Clin Gastroenterol Hepatol.* 2010;8(2):137-142.
3. de Paulo GA, Martins FP, Macedo EP, Goncalves ME, Mourao CA, Ferrari AP. *Endosc Int Open.* 2015;3(1):E7-E13.
4. Maitra S, Khanna P, Baidya DK. *Eur J Anaesthesiol.* 2014;31(5):266-273.
5. Eschertzhuber S, Brimacombe J, Hohlieder M, Keller C. *Anaesthesia.* 2009;64(1):79-83.<sup>1</sup>
6. Wong JG, Heaney M, Chambers NA, Erb TO, von Ungern-Sternberg BS. *Pediatr Anesth.* 2009;19(5):464-469.
7. van Zundert T, Brimacombe J.<sup>1</sup> *Anaesthesiol Intensive Ther.* 2012;44(2):63-70.<sup>1</sup>
8. Burgard G, Mollhoff T, Prien T. *J Clin Anesth.* 1996;8(3):198-201.
9. Chantzara G, Stroumpoulis K, Alexandrou N, Kokkinos L, Iacovidou N, Xanthos T. *Minerva Anesthesiol.* 2014;80(5):547-555.
10. Seet E, Yousaf F, Gupta S, Subramanyam R, Wong DT, Chung F. *Anesthesiology.* 2010;112(3):652-657.<sup>‡</sup>
11. Terblanche NCS, Middleton C, Choi-Lundberg DL, Skinner M. *Br J Anaesth.* 2018;120(2):353-360.

<sup>1</sup> Research sponsored by Teleflex Incorporated.

<sup>‡</sup> Research sponsored in part by Teleflex Incorporated.

Teleflex, the Teleflex logo, Cuff Pilot, LMA, and LMA Gastro are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. Revised: 02/2022. © 2022 Teleflex Incorporated. All rights reserved. MCI-2017-0748 · EN · REV 0 · 02 22 PDF

**Teleflex**

Distributed by: Teleflex Headquarters International, Ireland · Teleflex Medical Europe Ltd. · IDA Business & Technology Park  
Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · [orders.intl@teleflex.com](mailto:orders.intl@teleflex.com)  
United Kingdom Tel. +44 (0)14 94 53 27 61 · [info.uk@teleflex.com](mailto:info.uk@teleflex.com)  
South Africa Tel. +27 (0)11 807 4887 · [assist.africa@teleflex.com](mailto:assist.africa@teleflex.com)  
[www.teleflex.com](http://www.teleflex.com)